

ALT-R004 Three-Wheel Lightweight Walker with Seat & Bag



The **Alerta Three-Wheel Lightweight Walker** is equipped with a seat, backrest and storage pockets, and features a triangular foldable frame and adjustable handle height. It has robust 8" wheels and a lightweight frame for easy transport and manoeuvrability.

TECHNICAL SPECIFICATION

- > Handle height: 82-99cm
- > Seat height: 50cm
- > Max user weight: 135kg
- > Product weight: 6.25kg
- > Handle spacing: 33cm
- > Wheel size: 20cm (8")
- > Space between rear wheels: 58cm
- > Space between rear and front wheels: 49cm

KEY FEATURES & BENEFITS

- > Seat with backrest for users to rest
- > Robust lightweight frame
- > Lightweight and foldable for manoeuvrability and transport
- > Handle lock safety feature
- > Storage pockets
- > Robust 8" wheels for indoor and outdoor use
- > Reflectors for added visibility during low light