

ALT-R008 Alerta Aluminium Lightweight 2-in-1 Wheelchair Rollator



The **Alerta Aluminium Lightweight 2-in-1 Wheelchair Rollator** is an excellent dual purpose, hybrid mobility companion.

Its adaptability enables use as a sturdy lightweight rollator promoting independence, exercise and well-being, but also easily converts into a wheelchair by simply swinging the footrests into place and flipping the backrest over. It is the perfect friend promoting mobility by enabling the user to walk as far as they wish to without restricting distance from home, because it can then become a wheelchair and the journey can continue with the help of an assistant.

Alerta's Wheelchair Rollator surrounds the user with luxurious comforts such as soft extra-comfort padded seat, generously sized sturdy comfort backrest, quick release adjustable height ergonomic shaped handles, pull up break & push down hand-brake lock, easy swivel shock-absorbing wheels, storage bag with zip, swing-away detachable footrests, fixed padded armrests for wheelchair position, and easy one-click pull rapid fold for storage and transport.

KEY FEATURES & BENEFITS

- > Adjustable height handles
- > Adjustable back wheels
- > Spacious storage bag with detachable shoulder strap
- > Easy quick-fold design
- > Convert rollator to wheelchair
- > Walking stick holder
- > Seat with adjustable back rest
- > Detachable and adjustable footrests
- > Reflectors on front wheels
- > Handle lock safety feature
- > Robust wheels for indoor and outdoor use

TECHNICAL SPECIFICATION

- > Length (footrests stored): 75cm
- > Length (footrests deployed): 105cm
- > Width: 64cm
- > Seat height: 55cm
- > Seat width: 46cm
- > Adjustable handle height: 86-100cm (6 positions)
- > Max user weight: 136kg
- > Storage bag capacity: 5kg
- > Product weight: 9.9kg