

Four-Wheel Aluminium Lightweight Indoor Rollator



Blue **ALT-R010/B**



Green **ALT-R010/G**



Red **ALT-R010/R**

The **Alerta Four-Wheel Aluminium Lightweight Indoor Rollator** is an excellent personal home-help friend, equipped with a useful tray, bag, pull up break & push down hand-brake lock, easy swivel wheels, and features a robust, foldable frame and adjustable handle height. The lightweight narrow construction design offers easy transportation, and the sturdy aluminium frame gives stability and support when in use.

TECHNICAL SPECIFICATION

- > Handle height: 89-93.5cm
- > Max user weight: 120kg
- > Product weight: 5.24kg
- > Product dimensions: 57w x 65l
- > Tray size: 38 x 29cm
- > Basket size: 37w x 27l x 18d cm
- > Handle width (widest point): 41cm
- > Handle width (narrowest point): 28cm
- > Wheel size: 15cm
- > Space between front wheels: 44.5cm
- > Space between rear wheels: 48.5

KEY FEATURES & BENEFITS

- > Available in three colours (blue, red, green)
- > Robust lightweight frame
- > Easy-fold for storage and transport
- > Handle lock safety feature
- > Storage basket
- > Contoured transport tray
- > Adjustable height handles
- > Robust 6" wheels for indoor use