

ALT-R014 Three-Wheel Lightweight Walker with Basket, Tray & Bag



The **Alerta Three-Wheel Lightweight Walker** is equipped with a basket, tray and bag, and features a triangular foldable frame and adjustable handle height. It has robust 8" wheels and a lightweight frame for easy transport and manoeuvrability.

TECHNICAL SPECIFICATION

- > Handle height: 82-97cm
- > Max user weight: 115kg
- > Product weight: 6.5kg
- > Handle spacing: 45cm
- > Max width: 65cm
- > Wheel size: 20cm (8")
- > Tray/Basket size: 28 x 24cm
- > Basket height: 19cm
- > Max basket load: 5kg

KEY FEATURES & BENEFITS

- > Robust lightweight frame
- > Lightweight and foldable for manoeuvrability and transport
- > Handle lock safety feature
- > Spacious basket and storage bag
- > Includes tray with raised edges
- > Robust 8" wheels for indoor and outdoor use