

ALT-R002 Four-Wheeled Rollator with Solid Back



The **Alerta Four-Wheeled Aluminium Rollator with Solid Back** is a feature rich walking aid to assist users both around the house and when out and about.

It boasts a strong lightweight aluminium frame equipped with a padded seat for that essential rest, a bag and sturdy backrest, puncture proof tyres and robust build for the daily routines.

The ergonomic handgrips on easy to adjust height adjustable handles ensure it is perfectly suited to extended range of age groups with varying strength and abilities of hand use and grip power.

Alerta's lightweight Rollator has an easy folding mechanism for quick, convenient storage.

TECHNICAL SPECIFICATION

- > Handle height: 82 - 91cm
- > Seat height: 52cm
- > Max user weight: 136kg
- > Product weight: 6.5kg

KEY FEATURES & BENEFITS

- > Seat with sturdy backrest for users to rest
- > Robust lightweight aluminium frame
- > Quick fold design for easy transport and maneuverability
- > Handle lock safety feature
- > Spacious storage bag
- > Robust wheels for indoor and outdoor use
- > Detachable and foldable backrest
- > Hinged seat
- > Adjustable height handles