

BBS Revolution™



The easiest to use bladder scanner on the market

BBS Revolution™ Features

- 8 Transducer Probe (8 ultrasound signals) allows:-
- Automatic location of bladder without aiming.
- Automatic detection of male or female anatomy.
- Rapid (0.2 secs) and precise bladder volume measurement.
- Automatic self calibration - no annual service costs.
- Full 3 year warranty, includes accidental damage to probe.
- Includes Unlimited Certificated Remote Training.

BBS Revolution™ on Trolley




BBS™ Revolution on Desk Top Stand



BBS Revolution™ Probe




One Transducer position 1




Typical initial one transducer scanner position and estimate of resultant ultrasound picture.




One Transducer position 1




Screen may show off centre positioning.




One Transducer position 2




Once the ideal position is identified, scanner will provide bladder volume.




BBS Wide Angle Scan



BBS Revolution is the first scanner with a 180° field of view. This allows the scanner to automatically locate the bladder instantly. The sample ultrasound image shows the ease of locating the bladder.



BBS Wide Angle Scan



The BBS Revolution automates the scanning process, allowing users with very little experience and training to make a successful scan.

