

beurer

HealthManager

Modern health management. All the Beurer Connect products are remarkable because they are not only save and display measurements on the device itself; users can also transfer all the data to a PC or gain an overview on the move with a smartphone or mobile device.



The HealthManager represents the core system. All products carrying the HealthManager logo include the HealthManager software and an app version as a free download. There are different types of data transfer depending on the product type: Data is transmitted directly to your smartphone via NFC (near field communication), via *Bluetooth*® technology or to your PC using the conventional method, by USB. You can save the data locally or synchronise your data with any of your systems using the web version at any time – compatible with PC, iOS or Android.



Data Privacy/
Data Security
www.tuv.com
ID 000042517

Data security and privacy protection. TÜV Rheinland i-sec GmbH certified the Beurer Cloud „HealthManager“ according to „Data security and privacy protection of online applications“.



¹ Audited online application Beurer HealthManager cloud
² *Bluetooth*®: Data transfer with iOS and Android version of the HealthManager App
³ NFC: Data transfer with Android version of the HealthManager App



Beurer HealthManager Products overview



Connection via *Bluetooth*®



Connection via NFC



Connection via USB



756.636 - 0816 Subject to errors and changes

Beurer GmbH, P.O. Box 1427, 89004 Ulm, Germany, Soeflinger Strasse 218, 89077 Ulm, Germany
Tel. +49 731/39 89-0, Fax. +49 731/39 89-139, export@beurer.de www.beurer.com www.beurer-healthguide.com



beurer

HealthManager

Innovative technologies for your health.
The new generation of linked product world.



Data Privacy/
Data Security
www.tuv.com
ID 000042517




¹ Audited online application Beurer HealthManager cloud



Health Monitoring System

Retrieve personal data and evaluation **safe, fast, any time and anywhere**



Data available via web and app.
Access it on the move – any time and anywhere.



¹ Audited online application Beurer HealthManager cloud

HealthManager app

Connection via *Bluetooth®* or NFC technology



HealthManager PC software

The free **HealthManager PC software** enables you to analyse your body values even more closely. Additional views and settings allow for individual evaluations. You can use export or email functions to conveniently export your data and share it with others. You can also easily synchronise your data within the HealthManager system – **ONE complete system to monitor your health.**

Connection via USB interface



The Beurer HealthManager – the core of the Beurer Connect.

- **Free** We will offer you the HealthManager system as a free download for all products indicated as HealthManager products
- **Easy to use** and easy to install user-friendly interface
- **Fast and convenient** The application is designed specifically to suit your needs – with a time-saving and convenient display
- **At home and on the move** Your personal evaluations and data can be accessed any time and anywhere



HealthManager web

You can always keep an eye on your data with the free **HealthManager web** version. The exchange of your health data between HealthManager app and HealthManager PC software will take place via the cloud. You will find the online portal HealthManager web at connect.beurer.com

The app and software you download may differ from the illustrations.